
# IFU-UGMO

**ENG** English

## Universal Gutter Multifunction Orthosis

### Instructions for Use

Before using the device, please read the following instructions completely and carefully. Correct application and care is vital to the proper functioning of the device.

#### INTENDED USE / INDICATIONS:

- Treatment of minimally-displaced fractures of the metacarpals, sprains and strains, tendonitis and other cumulative trauma injuries.
- Provide support in the post-healing rehab phase.

#### CONTRAINDICATIONS:

Complex, compound and/or comminuted fractures.

#### WARNINGS AND PRECAUTIONS:

- Device to be fitted by a qualified healthcare professional.
- Patients should consult their medical professional immediately if they experience any pain, swelling, sensation changes, or any unusual reactions while using this device.
- This Device should not be used on a patient with polyester or elastane allergy.
- Ensure the orthosis is not too hot prior to application.
- Do not apply over open wound.
- Instruct the patient on correct adjustment of straps to maintain snug fit.
- The typical flexed position of the fingers in a gutter splint may cause skin maceration of the palm and fingers if the cast remains wet for more than several hours. Advise the patient to allow the cast to dry thoroughly before getting it wet again.
- Remove the device before any re-molding with a heat-gun, to avoid burning the patient.
- The orthosis will lose its shape above 57°C/135°F or under extreme sunlight. Cover with stockinet or clothing.
- Avoid water or vapor above 40°C/104°F for prolonged periods. Avoid heat sources (e.g. ovens, flames, sunny car windows, very humid conditions) for long periods.
- This device reduces mobility & dexterity. Do not operate heavy equipment.
- Patients should consult their medical professional before traveling by airplane.

#### 1. SELECT & GAUGE (See matching photos on packaging insert)

- Select orthosis size based on sizing chart on package.
- Select radial or ulnar treatment configuration. Conversion from ulnar to radial configuration is achieved by pushing out the thumbhole plug **F** and pushing fabric liner through the pre-cut thumbhole slot.

c) To pre-size the orthosis, offer the unpacked device to the patient's arm on the injured side and align the distal edge of the device with the desired end point on the finger. Assess the degree of wrist coverage obtained and note what degree of shortening, if any, is required.

#### 2. HEAT

• **Convection oven:** Place the orthosis (fabric side down) into a standard convection oven at 75-85°C (167-185°F) for 9 minutes or until it becomes soft / malleable. DO NOT heat for a total time greater than one hour.

• **Fastform or Nu-Wave IR Oven (1500W):** Lay the orthosis on the FastForm rack, with the fabric down ensuring the plastic does not touch the sides. Heat UGMO orthosis at **60% Power** for **5 mins** or per the specific settings on your appliance. Check that all areas are fully malleable. If not fully malleable, re-heat in 1-minute increments.

• The UGMO may be heated by placing the **unopened** package (fabric side up) in a waterbath at 75-85degC for 3-7 minutes or until fully malleable. Note: to prevent burning the patient, do not submerge the orthosis in hot water unless placed in a watertight poly bag.

#### 3. FIT

- Trim, if required, at the trim sections with scissors while the orthosis is hot.
- Fold and stack the last level onto the next level and compress to form a comfortable cuff.

#### 4. APPLY

- Place required finger or fingers in the finger sock **D** and close the wrist section of the device using Easy-fix tab **A**.
- Close the finger section of the device using the second Easy-fix tab **B**. Ensure the velour section **E** of the finger closure system sits between the immobilised and mobile fingers.
- Confirm good alignment and positioning of the device prior to bandage wrap.

#### 5. MOLD

- Wrap an elastic bandage (not included) firmly but not tightly from distal to proximal end to conform the entire orthosis. A damp bandage accelerates cooling.
- Manipulate the device into the most appropriate position for healing.

#### 6. SET

- Rest patient's arm on a table until the orthosis is rigid then remove the elastic bandage.

#### 7. SECURE

- Apply the Compliance-Lock strap to the wrist area of the orthosis. Thread it through the D-ring & tighten for a snug fit.
- To lock this strap, push the white fastener peg through the nearest matching hole on the strap and snap to lock in place. Thus the strap can be loosened a pre-determined distance e.g. to allow for swelling. To tighten it further break link between fastener and D-ring. To remove the compliance strap, cut the white peg with scissors.

Note: the physician should advise the patient on how to adjust the straps for a snug fit. Note: After removal the orthosis can be reformed with a heat gun in a localised area or in an oven. See [www.FastFormMedical.com](http://www.FastFormMedical.com) for further details.

#### 8. REHAB

- To prepare orthosis for rehab use, remove it & flex both angled hinges **G** fully for easy opening.

#### GETTING IT WET

If allowed by the physician, it is acceptable for the orthosis to become wet during use. Note: Where an acute finger position is molded in this may trap moisture and thus limit the drying properties of the device. The patient should allow the cast to dry fully before getting it wet again. It may take a number of hours to become completely dry. If allowed by the physician to loosen or remove the orthosis, the patient may dry it using a hair dryer on a low-heat setting. The physician may therefore advise limitations to the wetting of the device

#### CLEANING INSTRUCTIONS:

To clean the orthosis, rinse it in warm water. Note: it is not recommended to use household cleaners or personal hygiene products directly on the orthosis.

#### STORAGE:

This device is not sterile and should be stored between 4 - 35°C (39-95°F). It must also be protected from strong direct sunlight.

**WARRANTY:** FastForm Research Ltd covers defects in workmanship and materials for up to 12 weeks after this device is dispensed. Damage caused by normal wear & tear, abuse or alterations is not covered.

For further details on heating & application method, see: [www.FastFormmedical.com](http://www.FastFormmedical.com)

- Custom mouldable for single patient use.
- Not made with natural rubber latex -Rx only

**Notice:** while every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this device

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**FASTFORM** FastForm Research Ltd.,
**BREG**® Breg, Inc.,