PS-2603-01 rev H1 Feb 2018

## Instructions for Use Short Arm Multifunction Orthosis

**Before using the device, please read the following instructions completely and carefully. Correct application and care is vital to proper functioning of the device.**

### INTENDED USE / INDICATIONS:

- Treatment of minimally-displaced fractures of the wrist, sprains and strains, tendonitis and other cumulative trauma injuries.
- Provide support in the post-healing rehab phase.

### CONTRAINDICATIONS:

Complex, compound and/or comminuted fractures.

### WARNINGS AND PRECAUTIONS:

Device to be fitted by a qualified healthcare professional.

Patients should consult their medical professional immediately if they experience any pain, swelling, sensation changes, or any unusual reactions while using this device.

This Device should not be used on a patient with polyester or elastane allergy.

Ensure the orthosis is not too hot prior to application.

Do not apply over open wound.

Instruct the patient on correct adjustment of straps to maintain snug fit.

Remove the device before any re-molding with a heat-gun, to avoid burning the patient.

The orthosis will lose its shape above 57°C/135°F or under extreme sunlight. Cover with socknet or clothing.

Avoid water or vapor above 40°C/104°F for prolonged periods. Avoid heat sources (e.g. ovens, flames, sunny car windows, very humid conditions) for long periods.

This device reduces mobility & dexterity. Do not operate heavy equipment.

Patients should consult their medical professional before traveling by airplane.

Not suitable for children under 3 years of age due to small parts & choking hazard.

Remove white buttons on compliance strap, if not used, to avoid choking hazard.

### 1. SELECT & GAUGE (See matching photos on packaging insert)

a) Select orthosis size based on sizing chart on package.

b) To pre-size the orthosis, place it on the patient's arm, aligning the distal dorsal end parallel to the knuckles and measure from the end of the orthosis to the crease on the inside of the elbow. It should be approximately 3 finger-widths away from the crease. Note how many strips (if any) need to be cut off after heating.

### 2. HEAT

**Convection oven:** Place the orthosis (fabric side down) into a standard convection oven at 75-85°C (167-185°F) for 9 minutes or until it becomes soft / malleable. DO NOT heat for a total time greater than one hour.

**Fastform or Nu-Wave IR Oven (1500W):** Lay the orthosis on the FastForm rack, with the fabric down and the thumb section flat, ensuring the plastic does not touch the sides. Heat SAMO orthosis at **60% Power for 5 mins** or per specific settings on your appliance. Check that all areas are fully malleable. If not fully malleable, re-heat in 1-minute increments.

**Waterbath:** The SAMO may be heated by placing the **unopened** package (fabric side up) in a waterbath at 75-85°C (167-185°F) for 3-7 minutes or until fully malleable. To prevent burning the patient, do not submerge the orthosis in hot water unless placed in a watertight poly bag.

### 3. FIT

- Trim, if required, at the trim sections **C** with scissors while the orthosis is hot.
- Fold and stack the last level onto the next level and compress to form a comfortable cuff.

### 4. APPLY

- Rest the patient's elbow on a flat surface with the thumb abducted. Fit the orthosis over the thumb and align the radial hinge **D** along the radius. Ensure the thumbhole bridge **I** is snugly positioned.
- Stretch & wrap the palmar cuff **E** around the hand and align the closure tab **G** onto the closure fixing zone **F** & press into place.
- Stretch & wrap wing **B** around the arm to overlap wing **A** over its full length. Align & secure yellow dot **N** with yellow dot **M**. Ensure semi-circular tab **H** overlaps palmar area **J**.
- Easy-fix tabs **L** will temporarily secure wing in place. Confirm the radial hinge **D** is still aligned along the radius.

### 5. MOLD

- Wrap an elastic bandage (not included) firmly but not tightly to conform the entire orthosis. A bandage wetted with iced water accelerates cooling.
- Manipulate the device into the most appropriate position for healing, typically using Arm-wrestler grip and cupping for dorsal support up the forearm proximally.

### 6. SET

- Rest patient's arm on a table until the orthosis is rigid. Then remove the elastic bandage and the Easy-fix tabs. Verify that the medial wing **B** and palmar cuff **E** are free to allow for subsequent adjustment.

### 7. SECURE

- Apply the Compliance-Lock strap to the wrist area of the orthosis. Thread the strap through the slot in the liner. Thread it through the D-ring & tighten for a snug fit.
- To lock this strap, push the white fastener peg through the nearest matching hole on the strap and snap to lock in place. Thus the strap can be loosened a pre-determined distance e.g. to allow for swelling. To tighten it further break link between fastener and D-ring. To remove the compliance strap, cut the white peg with scissors.
- Fit forearm tension strap and tension for a snug fit.

Note: the physician should advise the patient on how to adjust the straps for a snug fit. Note: After removal the orthosis can be reformed with a heat gun in a localised area or in an oven. See [www.FastFormMedical.com](http://www.FastFormMedical.com) for further details.

### 8. REHAB

- To prepare orthosis for rehab use, remove it & flex radial hinge **D** fully for easy opening.
- Likewise flex the Palmar Cuff Hinge.

### GETTING IT WET

If allowed by the physician, it is acceptable for the orthosis to become wet during use. It may take a number of hours to become completely dry. If allowed by the physician to loosen or remove the orthosis, the patient may dry it using a hair dryer on a low-heat setting.

### CLEANING INSTRUCTIONS:

To clean the orthosis, rinse it in warm water. Note: it is not recommended to use household cleaners or personal hygiene products directly on the orthosis.

### STORAGE:

This device is not sterile and should be stored between 4 - 35°C (39-95°F). It must also be protected from strong direct sunlight.

**WARRANTY:** FastForm Research Ltd covers defects in workmanship and materials for up to 12 weeks after this device is dispensed. Damage caused by normal wear & tear, abuse or alterations is not covered.

For further details on heating & application method, see:

[www.breg.com/fastform/HeatingTimes](http://www.breg.com/fastform/HeatingTimes) or [www.FastFormmedical.com](http://www.FastFormmedical.com)

**-Custom mouldable for single patient use.**

**-Not made with natural rubber latex -Rx only**

**Notice:** while every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this device.

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**Breg, Inc.,** 2885 Loker Ave. East, Carlsbad, CA 92010, USA