**ENG VERSI™ Sheet Instructions for Use**

Before using the device, please read the following instructions completely and carefully. Correct application and care is vital to the proper functioning of device.

**INTENDED USE / INDICATIONS:**

The sheet is indicated for use by a trained healthcare professional in the custom crafting of bespoke orthopedic supports including extensions to other FastForm® Products; hand, wrist, finger and thumb splints; orthopedic orthosis, external immobilization devices and rehabilitation aids; serial splints; outrigger bases; spasticity splints; functional position splints, resting mitt splints.

**CONTRAINDICATIONS:**

Complex, compound and/or comminuted fractures and load bearing applications.

**WARNINGS AND PRECAUTIONS:**

- The VERSI™ Sheet is designed for use by a qualified healthcare professional.
- Patients should consult a medical professional immediately if they experience any pain, swelling, sensation changes, or unusual reactions while using their device.
- This Device should not be used on a patient with polyester or elastane allergy.
- Ensure the crafted device is not too hot prior to application.
- Do not apply over open wound.
- Remove the device before re-molding with heat-gun, to avoid burning patient.
- The crafted device will lose its shape above 57°C/135°F or under extreme sunlight. Cover with stockinet or clothing.
- Avoid water or vapor above 40°C/104°F for prolonged periods. Avoid heat sources (e.g. ovens, flames, sunny car windows, very humid conditions) for long periods.

**HEAT**

- **Convection oven:** Place the VERSI™ Sheet (fabric side down) into a standard convection oven at 75-85°C (167-185°F) for 9 minutes or until it becomes soft / malleable. DO NOT heat for a total time greater than one hour.
- **FastForm® or Nu-Wave IR Oven (1500W):** Lay the VERSI™ Sheet on the FastForm® rack, with the fabric down, ensuring the plastic does not touch the sides. Heat VERSI™ Sheet at 60% Power for 5 mins. Check that all areas are fully malleable.
- If not fully malleable, re-heat in 1-minute increments.
- The VERSI™ Sheet may be heated by placing the **unopened** package (fabric side up) in a water bath at 75-85degC for 3-7 minutes or until fully malleable. Note: to prevent burning the patient, do not submerge the sheet in hot water unless placed in a watertight poly bag.
- Remove backing paper from sheet. These areas can be used for tessellation of multiple sheets.

**FIT:** a) Trim the sheet to the desired shape, with scissors while the sheet is hot. b) Fold, stack and compress to form a comfortable cuff.

**APPLY:** Ensure that the fabric side of the sheet is applied to the patient's skin. The plastic side should never touch the skin.

**GETTING IT WET**

If allowed by the physician, it is acceptable for the device to become wet during use. It may take a number of hours to become completely dry. If allowed by the physician to loosen or remove the crafted device, the patient may dry it using a hair dryer on low-heat setting. To avoid maceration it must dry fully before getting wet again.

**CLEANING INSTRUCTIONS:**

To clean the crafted device, rinse it in warm water. Note: it is not recommended to use household cleaners or personal hygiene products directly on the crafted device.

**STORAGE:** This VERSI™ Sheet is not sterile and should be stored between 4 - 35°C (39-95°F). It must also be protected from strong direct sunlight.

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
FastForm Research Ltd.  
Unit 4, Lear Avenue,  
Airport Business Park,  
Killowen,  
Co. Waterford,  
Ireland  
[www.fastformmedical.com](http://www.fastformmedical.com)


Not made with natural rubber latex

**WARRANTY:** FastForm Research Ltd covers defects in workmanship and materials for up to 12 weeks after this VERSI™ Sheet is dispensed. Damage caused by normal wear & tear, abuse or alterations is not covered.

**Notice:** While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

82% Polytrex® Polymer, 14% Polyester, 3% Polyamide, 1% Elastane  
- Not made with natural rubber latex

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