ENG VERSI™ Sheet Instructions for Use

Before using the device, please read the following instructions completely and carefully. Correct application and care is vital to the proper functioning of device. **INTENDED USE / INDICATIONS:**

The sheet is indicated for use by a trained healthcare professional in the custom crafting of bespoke orthopedic supports including extensions to other FastForm® Products; hand, wrist, finger and thumb splints; orthopedic orthosis, external immobilization devices and rehabilitation aids; serial splints; outrigger bases; spasticity splints; functional position splints, resting mitt splints.

CONTRAINDICATIONS:

Complex, compound and/or comminuted fractures and load bearing applications.

WARNINGS AND PRECAUTIONS:

- The VERSI™ Sheet is designed for use by a qualified healthcare professional.
- Patients should consult a medical professional immediately if they experience any pain, swelling, sensation changes, or unusual reactions while using their device.
- This Device should not be used on a patient with polyester or elastane allergy.
- Ensure the crafted device is not too hot prior to application.
- Do not apply over open wound.
- Remove the device before re-molding with heat-gun, to avoid burning patient.
- The crafted device will lose its shape above 57°C/135°F or under extreme sunlight. Cover with stockinet or clothing.
- Avoid water or vapor above 40°C/104°F for prolonged periods. Avoid heat sources (e.g. ovens, flames, sunny car windows, very humid conditions) for long periods.

HEAT

- Place the VERSI™ Sheet (fabric side down) in the heater. Ensure the plastic does not touch the sides. Check that all areas are fully malleable. If not, reheat in 1 minute increments. DO NOT heat for a total time greater than one hour.

- **Convection oven:** = 9 mins in pre-heated oven, 75-85°C (167-185°F)

- **Nu-Wave Pro Plus / Elite:** = 5 mins @ 220°F

- **Nu-Wave Primo:** = 5 mins @ 260°F

- **Water bath:** Place the unopened package (fabric side up) in a water bath at 75-85°C (167-185°F) for 3-7 minutes or until fully malleable. To prevent burning the patient, do not submerge the sheet in hot water unless placed in a watertight poly bag.

- Remove backing paper from sheet. These areas can be used for tessellation of multiple sheets.

FIT: a) Trim the sheet to the desired shape, with scissors while the sheet is hot.

b) Fold, stack and compress to form a comfortable cuff.

APPLY: Ensure that the fabric side of the sheet is applied to the patient's skin. The plastic side should never touch the skin.

GETTING IT WET

If allowed by the physician, it is acceptable for the device to become wet during use. It may take a number of hours to become completely dry. If allowed by the physician to loosen or remove the crafted device, the patient may dry it using a hair dryer on low-heat setting. To avoid maceration it must dry fully before getting wet again.

CLEANING INSTRUCTIONS:

To clean the crafted device, rinse it in warm water. Note: it is not recommended to use household cleaners or personal hygiene products directly on the crafted device.

STORAGE: This VERSI™ Sheet is not sterile and should be stored between 4 - 35°C (39-95°F). It must also be protected from strong direct sunlight.

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
FastForm Research Ltd.
Unit 4, Lear Avenue,
Airport Business Park,
Killowen,
Co. Waterford,
Ireland
www.fastformmedical.com


Not made with natural rubber latex

WARRANTY: FastForm Research Ltd covers defects in workmanship and materials for up to 12 weeks after this VERSI™ Sheet is dispensed. Damage caused by normal wear & tear, abuse or alterations is not covered.

Notice: While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product. If a serious incident occurs with this device, please report it to the Manufacturer and your Competent Authority.

82% Polytrex® Polymer, 14% Polyester, 3% Polyamide, 1% Elastane
- Not made with natural rubber latex

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